



Dimensions (inches)

Ref.	External dimensions			Internal dimensions - without work surface		
	Width	Depth	Height*	Width	Depth	Height
Sd 48 I	50 3/16"	28 7/16"	45 2/8" min 52 5/8" max	48 13/16"	23 1/4"	35"

*According to filtration column type

Specifications

Number of columns	1
Number of fans (IP44)	1
Processed air flow	135 cfm
Air velocity at the openings in working position	79 fpm to 118 fpm
Voltage / Frequency	85-260 V / 50-60 Hz
Total power consumption	45 watts (without lighting)
Amperage absorbed	1,6 amp.
Sound level	52 dbA