



Dimensions (inches)

Ref.	External dimensions			Internal dimensions - without work surface		
	Width	Depth	Height*	Width	Depth	Height
XL 483	50 3/16"	37"	45 2/8" min 52 5/8" max	48 13/16"	31 1/4"	35"

*According to filtration column type

Specifications

Number of columns	3
Number of fans (IP44)	3
Processed air flow	405 cfm
Air velocity at the openings in working position	79 fpm to 118 fpm
Voltage / Frequency	85-260 V / 50-60 Hz
Total power consumption	135 watts (without lighting)
Amperage absorbed	4,8 amp.
Sound level	58 dbA